

A proven non invasive technology

Diode laser uses infra red light to accelerate the production and remodeling of mucosa collagen.

Age and child birth and stress are most common factors responsible for muscular atrophy process

which results in dryness, sexual difficulties, itching, burning and painful Intercourse and tissue laxity as well urinary incontinence.

The main cause is loss of tone of vaginal mucosa which is targeted while doing treatment with diode laser which has bio stimulatory effect that stimulated new collagen formation and regeneration of epithelium.

## Laser Vaginal Rejuvenation

Improves muscle tone

Highly effective & Painless

Non- invasive

Enhanced sexual gratification

Effective in vaginal tightening

Improved sensation

Increased lubrication and moisture

## Symptom

Vaginal laxity

Stress Urinary Incontinence

Vaginal Dryness

Vaginal Infection

Post Delivery Rehabilitation

Polyp, Skin Tag Removal



Bringing Innovation Closer To You

Recreate

Rejuvenation

Revitalise



# TECHNICAL SPECIFICATIONS

| Specifications         | Aqua Plus           | Trident        |
|------------------------|---------------------|----------------|
| Maximum Power          | 15W + 200mW         | 15W+15W        |
| Wavelength             | 1470nm+650nm        | (980nm+1470nm) |
| Magic Pulse Technology | No                  | Yes            |
| Mode                   | Continuous & Pulsed |                |
| Aiming beam wavelength | 650nm               |                |
| Interlock Key          | Provided            |                |
| Optical Fiber          | Bare                |                |
| Fiber thickness        | 600 Microns         |                |
| Fiber Connector        | SMA 905A            |                |
| Fiber Length           | 3 Meters            |                |
| Display screen size    | 10.5 inches         |                |
| Foot switch            | Bluetooth & Wired   |                |
| Stripper & Cutter      | 1                   |                |
| Surgical Hand piece    | 1                   |                |

## Why Patient Prefer Diode Laser Treatment

Quick , Easy & Effective

Outpatient Treatment


No Cut

No stitches

No Anesthesia required

Absolutely Painless

Highly Effective



ISO 13485:2016



ISO 9001:2015

Marketed by  
**AXON Medical Solutions Pvt.Ltd.**

Unit No.11, Pace City – 1, Sector – 37, Honda Chowk, Gurgaon, Haryana - 122001

+91 9599017699 / 9810333158

info@axonmedicalsolutions.com

www.axonmedicalsolutions.com

Manufactured by  
**Intense Medical & Dental System Pvt. Ltd.**

G-7, Industrial Area Bagroda Bhopal M.P.- 462026

+91 9350014462

imdsl@yahoo.com

imdsl.co.in